



KACTUS

Cynomolgus/Rhesus macaque TNFSF15 Protein (FSF-CM115) 100UG

特徴

Description: Recombinant Cynomolgus/Rhesus macaque TNFSF15 Protein is expressed from Expi293 with His tag at the N-terminal. It contains Leu72-Leu251.



シリーズ共通仕様

種由来	Cynomolgus/Rhesus macaque	発現宿主	Expi293
純度	> 95% as determined by Tris-Bis PAGE	タグ	N-His
分子量	The protein has a predicted MW of 21.52 kDa. Due to glycosylation, the protein migrates to 25-35 kDa based on Tris-Bis PAGE result.	ターゲット	TNFSF15
Gene/Protein ID	A0A2K5UA22	エンドトキシン	Less than 1EU per µg by the LAL method.

配列	Leu72-Leu251	データ シート	https://static.kactusbio.com/Attachment/product/FSF-CM115/Liquid-FS-F-CM115.pdf
----	--------------	------------	---

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	83-1237-07	Cynomolgus/Rhesus monkey TNFSF15 Protein (FSF-CM115) 100 UG	FSF-CM115-100UG	1個	84,000円	
	83-1237-08	Cynomolgus/Rhesus monkey TNFSF15 Protein (FSF-CM115) 500 UG	FSF-CM115-500UG	1個	281,000円	