



KACTUS

Cynomolgus CD2/SRBC Protein (CD2-CM202) 100UG

特徴

Description: Recombinant Cynomolgus CD2/SRBC Protein is expressed from Expi293 with His tag at the C-terminal. It contains Lys25-Asp209.



シリーズ共通仕様

種由来	Cynomolgus	発現宿主	Expi293
純度	> 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC	タグ	C-His
分子量	The protein has a predicted MW of 22.49 kDa. Due to glycosylation, the protein migrates to 35-45 kDa based on Tris-Bis PAGE result.	ターゲット	CD2/SRBC
Gene/Protein ID	Q6SZ61	エンドトキシン	Less than 1EU per µg by the LAL method.

配列	Lys25-Asp209	データ シート	https://static.kactusbio.com/Attachment/product/CD2-CM202/Liquid-CD2-CM202.pdf
----	--------------	------------	--

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	83-1237-41	Cynomolgus CD2/SR BC Protein (CD2-CM2 02) 100UG	CD2-CM202-100UG	1個	84,000円	
	83-1237-42	Cynomolgus CD2/SR BC Protein (CD2-CM2 02) 500UG	CD2-CM202-500UG	1個	281,000 円	