



KACTUS

Biotinylated Mouse GUCY2C/Guanylyl cyclase C Protein (GCC-MM401B) 100UG

特徴

Description: Recombinant Biotinylated Mouse GUCY2C /Guanylyl cyclase C Protein is expressed from Expi293 with His tag and Avi tag at the C-terminal. It contains Val20-Met433.



シリーズ共通仕様

種由来	Mouse	発現宿主	Expi293
純度	> 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC	タグ	C-His-Avi
分子量	The protein has a predicted MW of 50.01 kDa. Due to glycosylation, the protein migrates to 70-100 kDa based on Tris-Bis PAGE result.	ターゲット	GUCY2C/Guanylyl cyclase C
Gene/Protein ID	Q3UWA6-1	エンドトキシン	Less than 1EU per µg by the LAL method.

配列	Val20-Met433	データ シート	https://static.kactusbio.com/Attachment/product/GCC-MM401B/Liquid-GCC-MM401B.pdf
----	--------------	------------	---

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	83-1237-45	Biotinylated Mouse GUCY2C/Guanylyl c yclase C Protein (G CC-MM401B) 100U G	GCC-MM401B-100UG	1個	141,000 円	
	83-1237-46	Biotinylated Mouse GUCY2C/Guanylyl c yclase C Protein (G CC-MM401B) 500U G	GCC-MM401B-500UG	1個	467,000 円	