



KACTUS

Cynomolgus LY75/CD205 Protein (LY7-CM105) 100UG

特徴

Description: Recombinant Cynomolgus LY75/CD205 Protein is expressed from Expi293 with His tag at the C-terminal. It contains Ser28-Asp1666.



シリーズ共通仕様

種由来	Cynomolgus	発現宿主	Expi293
純度	> 95% as determined by Tris-Bis PAGE; > 90% as determined by HPLC	タグ	C-His
分子量	The protein has a predicted MW of 190.17 kDa. Due to glycosylation, the protein migrates to 230-245 kDa based on Tris-Bis PAGE result.	ターゲット	LY75/CD205
Gene/ Protein ID	XP_005573337.2	エンドトキシン	Less than 1EU per µg by the LAL method.

配列	Ser28-Asp1666	データ シート	https://static.kactusbio.com/Attachment/product/LY7-CM105/Liquid-LY7-CM105.pdf
----	---------------	------------	--

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	83-1237-69	Cynomolgus LY75/CD2 05 Protein (LY7-CM1 05) 100UG	LY7-CM105-100UG	1個	84,000円	
	83-1237-70	Cynomolgus LY75/CD2 05 Protein (LY7-CM1 05) 500UG	LY7-CM105-500UG	1個	281,000 円	