



KACTUS

Human HSP70 Protein (HSP-HE170) 100UG

特徴

Description: Recombinant Human HSP70 Protein is expressed from E.coli with His tag at the N-terminal. It contains Ala2-Asp641 (Glu110Asp).



シリーズ共通仕様

種由来	Human	発現宿主	E.coli
純度	> 95% as determined by Tris-Bis PAGE;> 95% as determined by HPLC	タグ	N-His
分子量	The protein has a predicted MW of 72.2 kDa same as Tris-Bis PAGE result.	ターゲット	HSP70
Gene/ Protein ID	PODMV8-1	エンドトキシン	Less than 1EU per µg by the LAL method.

配列	Ala2-Asp641(Glu110 Asp)	データシート	https://static.kactusbio.com/Attachment/product/HSP-HE170/Liquid-HSP-HE170.pdf
----	-------------------------	--------	---

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
--------	--------	-----	----	----	------------------------------	----

	83-1237-71	Human HSP70 Protein (HSP-HE170) 100UG	HSP-HE170-100UG	1個	73,000円	
---	----------------------------	---------------------------------------	---------------------------------	----	---------	--

	83-1237-72	Human HSP70 Protein (HSP-HE170) 500UG	HSP-HE170-500UG	1個	241,000円	
---	----------------------------	---------------------------------------	---------------------------------	----	----------	--