



KACTUS

Human FLT3 Ligand Protein (FLT-HE03L) 500UG

特徴

Description: Recombinant Human FLT3 Ligand Protein is expressed from E.coli without tag. It contains Thr27-Ala 181.



シリーズ共通仕様

種由来	Human	発 現 宿主	E.coli
純度	> 95% as determined by Tris-Bis PAGE;> 95% as determined by HPLC	タグ	No Tag
分子量	The protein has a predicted MW of 17.61 kDa same as Tris-Bis PAGE result.	ターゲ ット	FLT3 Ligand
Gene/ Prote inID	P49771-1	エンド トキシ ン	Less than 1EU per µg by the LAL method.

配列	Thr27-Ala181	データ シート	https://static.kactusbio.com/Attachment/product/FLT-HE03L/Lyophilized-FLT-HE03L.pdf
----	--------------	------------	---

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
--------	--------	-----	----	----	------------------------------	----

	83-1237-75	Human FLT3 Ligand Protein (FLT-HE03L) 100UG	FLT-HE03L-100UG	1個	73,000円	
---	----------------------------	---	---------------------------------	----	---------	--

	83-1237-76	Human FLT3 Ligand Protein (FLT-HE03L) 500UG	FLT-HE03L-500UG	1個	241,000円	
---	----------------------------	---	---------------------------------	----	----------	--