



KACTUS

Human APOE3/Apolipoprotein E Protein (APO-HM203) 100UG

特徴

Description: Recombinant Human APOE3/Apolipoprotein E Protein is expressed from Expi293 with hFc tag at the N-terminal. It contains Lys19-His317.



シリーズ共通仕様

種由来	Human	発現宿主	Expi293
純度	> 95% as determined by Tris-Bis PAGE	タグ	N-hFc
分子量	The protein has a predicted MW of 59.76 kDa. Due to glycosylation, the protein migrates to 60-70 kDa based on Tris-Bis PAGE result.	ターゲット	APOE3/Apolipoprotein E
Gene/Protein ID	P02649	エンドトキシン	Less than 1EU per µg by the LAL method.

配列	Lys19-His317	データ シート	https://static.kactusbio.com/Attachment/product/APO-HM203/Liquid-APO-HM203.pdf
----	--------------	------------	--

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
--------	--------	-----	----	----	------------------------------	----

	83-1237-83	Human APOE3/Apolipoprotein E Protein (APO-HM203) 100UG	APO-HM203-100UG	1個	73,000円	
---	------------	--	-----------------	----	---------	--

	83-1237-84	Human APOE3/Apolipoprotein E Protein (APO-HM203) 500UG	APO-HM203-500UG	1個	241,000円	
---	------------	--	-----------------	----	----------	--