



KACTUS

Mouse Complement factor I Protein (CFI-MM101) 100UG

特徴

Description: Recombinant Mouse Complement factor I Protein is expressed from Expi293 with His tag at the C-terminal. It contains Arg19-Val603.



シリーズ共通仕様

種 由来	Mouse	発 現 宿主	Expi293
純度	> 95% as determined by Tris-Bis PAGE; > 95% as determined by HPLC	タグ	C-His
分 子 量	The protein has a predicted MW of 66.39 kDa. Due to autocatalytic cleavage and glycosylation, the protein migrates to 80-100 kDa on Tris-Bis PAGE result.	ターゲ ット	Complement factor I
Gene/ProteinID	Q61129	エンド トキ シン	Less than 1 EU per µg by the LAL method.

配列	Arg19-Val603	データシート	https://static.kactusbio.com/Attachment/product/CFI-MM101/Liquid-CFI-MM101.pdf
----	--------------	--------	---

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	83-1237-89	Mouse Complement factor I Protein (CFI-M101) 100UG	CFI-MM101-100UG	1個	84,000円	
	83-1237-90	Mouse Complement factor I Protein (CFI-M101) 500UG	CFI-MM101-500UG	1個	281,000円	