



KACTUS

Cynomolgus NTS1 Protein (NTS-CM101) 100UG

特徴

Description: Recombinant Cynomolgus NTS1 Protein is expressed from Expi293 with His tag at the C-terminal. It contains Ser24-Leu163.



シリーズ共通仕様

種由来	Cynomolgus	発 現 宿主	Expi293
純度	> 95% as determined by Tris-Bis PAGE	タグ	C-His
分子量	The protein has a predicted MW of 18.84 kDa. Due to glycosylation, the protein migrates to 20-25 kDa based on Tris-Bis PAGE result.	ターゲット	NTS1
Gene/ Protein ID	Q2PG57	エンドトキシン	Less than 1EU per μ g by the LAL method.

配列	Ser24-Leu163	データ シート	https://static.kactusbio.com/Attachment/product/NTS-CM101/Liquid-NTS-CM101.pdf
----	--------------	------------	---

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
--------	--------	-----	----	----	------------------------------	----

	83-1237-91	Cynomolgus NTS1 Protein (NTS-CM101) 10 OUG	NTS-CM101-100UG	1個	75,000円	
---	----------------------------	--	---------------------------------	----	---------	--

	83-1237-92	Cynomolgus NTS1 Protein (NTS-CM101) 50 OUG	NTS-CM101-500UG	1個	252,000円	
---	----------------------------	--	---------------------------------	----	----------	--