



KACTUS

Biotinylated Human KRAS G12A (HLA-A*11:01) Protein (MHC-HM441B) 100UG

特徴

Description: Recombinant Biotinylated Human KRAS G12A (HLA-A*11:01) Protein is expressed from Expi293 with His tag and Avi tag at the C-terminal. It contains Gly25-Thr305(HLA-A*11:01), Ile21-Met119(B2M) and VVVGAGVGK peptide.



シリーズ共通仕様

種由来	Human	発現宿主	Expi293
純度	> 95% as determined by Tris-Bis PAGE;> 95% as determined by HPLC	タグ	C-His-Avi
分子量	The protein has a predicted MW of 50.30 kDa. Due to glycosylation, the protein migrates to 52-62 kDa based on Tris-Bis PAGE result.	ターゲット	KRAS G12A (HLA-A*11:01)
Gene/Protein ID	AAV53343.1(HLA-A*11:01)&P61769(B2M)&VVVGAGVGK	エンドキソン	Less than 1EU per μg by the LAL method

配列	Gly25-Thr305(HLA-A*11:01), Ile21-Met119(B2M)and VVGAAGVGK peptide	データシート	https://static.kactusbio.com/Attachment/product/MHC-HM441B/Lyophilized-MHC-HM441B.pdf
----	---	--------	---

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	83-1237-95	Biotinylated Human KRAS G12A (HLA-A* 11:01) Protein (M HC-HM441B) 100U G	MHC-HM441B-100UG	1個	83,000円	
	83-1237-96	Biotinylated Human KRAS G12A (HLA-A* 11:01) Protein (M HC-HM441B) 500U G	MHC-HM441B-500UG	1個	334,000 円	