

US Biological (antibodypab)

84-1663-07 [取扱停止] TAT, ID (TAT, Tyrosine amino transferase, L-tyrosine:2-oxoglutarate aminotransferase) (APC) 200ul 042618-APC



アズワン品番 84-1663-07

型番 042618-APC

入り数 1個

標準価格 -円

取扱停止



特徴

TAT is a mitochondrial protein tyrosine aminotransferase which is present in the liver and catalyzes the conversion of L-tyrosine into p-hydroxyphenylpyruvate.

Applications: Suitable for use in Western Blot, Immunohistochemistry, Flow Cytometry (Not Tested), FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500 Immunohistochemistry: 1:50-100

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size: 200ul

Host: rabbit

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity: Human

Isotype: IgG

Calc Applications Abbrev: FLISA IHC WB

Calc Crossreactivity: Hu

Immunogen: TAT antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 285~315 amino acids from the Central region of human TAT.

Shelf Life: 6months

Swiss Prot Number: P17735

EU Commodity Code: 30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp