

US Biological (antibodypab)

84-1663-09 [取扱停止] NCOA7, NT (NCOA7, ERAP 140, ESNA1, Nuclear receptor coactivator 7, 140kD estrogen receptor-associated protein, Estrogen nuclear receptor coactivator 1) (APC) 200ul 038837-APC



アズワン品番 84-1663-09

型番 038837-APC

入り数 1個

標準価格 -円

取扱停止



特徴

- NCOA7 enhances the transcriptional activities of several nuclear receptors.
- It is involved in the coactivation of different nuclear receptors, such as ESR1, THRB, PPARG and RARA.
- Applications: Suitable for use in Western Blot, Immunohistochemistry, Flow Cytometry (Not Tested), FLISA
- Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500 Immunohistochemistry: 1:50-100
- Storage and Stability: Store product at 4°C in the dark.
- DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.
- Dilute required amount only prior to immediate use.
- Further dilutions can be made in assay buffer.
- Caution: APC conjugates are sensitive to light.
- For maximum recovery of product, centrifuge the original vial prior to removing the cap.
- Note: Applications are based on unconjugated antibody.

仕様

- Size: 200ul
- Host: rabbit
- Source Antibody: human
- Grade: Affinity Purified
- Purity: Purified by Protein A affinity chromatography.
- Form: Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).
- Specificity: Human

- Isotype: IgG
- Calc Applications Abbrev: FLISA IHC WB
- Calc Crossreactivity: Hu
- Immunogen: NCOA7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 82~112 amino acids from the N-terminal region of human NCOA7.
- Shelf Life: 6months
- Swiss Prot Number: Q8NI08
- EU Commodity Code: 30021010
- この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)
- 製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

- 【試薬に関するお問合せ】
- アズワン株式会社 試薬・プロセス材料グループ
- TEL: 06-6447-8641
- FAX: 06-6447-8642
- E-mail: Reagents@so.as-1.co.jp