

US Biological (antibodypab)

84-1663-20 [取扱停止] SPRR1A, CT (SPRR1A, Cornif
in-A, 19kD pancornulin, SPRK, Small proline-rich prot
ein IA) (APC) 200ul 042265-APC



アズワン品番 84-1663-20

型番 042265-APC

入り数 1個

標準価格 -円

取扱停止



特徴

Cross-linked envelope protein of keratinocytes.

It is a keratinocyte protein that first appears in the cell cytosol, but ultimately becomes cross-linked to membrane proteins by transglutaminase.

All that results in the formation of an insoluble envelope beneath the plasma membrane.

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu

Immunogen:SPRR1A antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 60~89 amino acids from the C-terminal region of human SPRR1A.

Shelf Life:6months

Swiss Prot Number:P35321

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp