



US Biological (antibodypab)

84-1663-21 [取扱停止] KIF6, NT (KIF6, C6orf102, Kinesin-like protein KIF6) (APC) 200ul 037449-APC



アズワン品番 84-1663-21

型番 037449-APC

入り数 1個

標準価格 -円

取扱停止



特徴

KIF6 belongs to the kinesin-like protein family.

There are four named isoforms.

The KIF6 719Arg variant has been associated with an increased risk of coronary disease and an increased response to statin therapy.

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human, mouse

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu Mo

Immunogen:KIF6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 41~70 amino acids from the N-terminal region of human KIF 6.

Shelf Life:6months

Swiss Prot Number:Q6ZMV9

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp