

US Biological (antibodypab)

84-1663-22 [取扱停止] RPL37, CT (RPL37, 60S ribosomal protein L37, G1.16) (APC) 200ul 041213-APC



アズワン品番 84-1663-22

型番 041213-APC

入り数 1個

標準価格 -円

取扱停止



特徴

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit.

Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins.

This gene encodes a ribosomal protein that is a component of the 60S subunit.

The protein belongs to the L37E family of ribosomal proteins.

It is located in the cytoplasm.

The protein contains a C2C2-type zinc finger-like motif.

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein G affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu

Immunogen:RPL37 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 67~96 amino acids from the C-terminal region of human RPL37 Antibody (C-term).

Shelf Life:6months

Swiss Prot Number:P61927

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp