



US Biological (antibodypab)

84-1663-23 [取扱停止] NUPL2, NT (NUPL2, CG1, Nucleoporin-like protein 2, NLP-1, NUP42 homolog, Nucleoporin hCG1) (APC) 200ul 039225-APC



アズワン品番	84-1663-23
型番	039225-APC
入り数	1個
標準価格	-円

取扱停止



特徴

Required for the export of mRNAs containing poly(A) tails from the nucleus into the cytoplasm.

In case of infection by HIV-1, it may participate in the docking of viral Vpr at the nuclear envelope.

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu

Immunogen:NUPL2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 94~124 amino acids from the N-terminal region of human NUPL2.

Shelf Life:6months

Swiss Prot Number:O15504

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp