

US Biological (antibodypab)

84-1663-34 [取扱停止] PEX12, ID (PEX12, PAF3, Peroxisome assembly protein 12, Peroxin-12, Peroxisome assembly factor 3) (APC) 200ul 039888-APC

**US Biological**  
Life Sciences

アズワン品番 84-1663-34

型番 039888-APC

入り数 1個

標準価格 -円

取扱停止

**US Biological**  
Life Sciences

## 特徴

- Peroxins (PEXs) are proteins that are essential for the assembly of functional peroxisomes.
- The peroxisome biogenesis disorders (PBDs) are a group of genetically heterogeneous autosomal recessive, lethal diseases characterized by multiple defects in peroxisome function.
- The peroxisomal biogenesis disorders are a heterogeneous group with at least 14 complementation groups and with more than 1 phenotype being observed in cases falling into particular complementation groups.
- Although the clinical features of PBD patients vary, cells from all PBD patients exhibit a defect in the import of one or more classes of peroxisomal matrix proteins into the organelle.
- Applications: Suitable for use in Western Blot, FLISA
- Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500
- Storage and Stability: Store product at 4°C in the dark.
- DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.
- Dilute required amount only prior to immediate use.
- Further dilutions can be made in assay buffer.
- Caution: APC conjugates are sensitive to light.
- For maximum recovery of product, centrifuge the original vial prior to removing the cap.
- Note: Applications are based on unconjugated antibody.

## 仕様

- Size: 200ul
- Host: rabbit
- Source Antibody: human
- Grade: Affinity Purified
- Purity: Purified by Protein A affinity chromatography.

- Form: Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).
- Specificity: Human
- Isotype: IgG
- Calc Applications Abbrev: FLISA WB
- Calc Crossreactivity: Hu
- Immunogen: PEX12 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 137~167 amino acids from the Central region of human PEX12.
- Shelf Life: 6months
- Swiss Prot Number: O00623
- EU Commodity Code: 30021010
- この商品は法規制を確認していません。(法規制によって販売できない場合もございます)
- 製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。
  
- 【試薬に関するお問合せ】
- アズワン株式会社 試薬・プロセス材料グループ
- TEL: 06-6447-8641
- FAX: 06-6447-8642
- E-mail: Reagents@so.as-1.co.jp