

US Biological (antibodypab)

84-1663-36 [取扱停止] ACSS2, NT (ACSS2, ACAS2, Acetyl-coenzyme A synthetase, cytoplasmic, Acetate-CoA ligase, Acetyl-CoA synthetase, Acyl-CoA synthetase short-chain family member 2, Acyl-activating enzyme) (APC) 200ul 031566-APC



アズワン品番	84-1663-36
型番	031566-APC
入り数	1個
標準価格	-円
取扱停止	



特徴

ACSS2 is a cytosolic enzyme that catalyzes the activation of acetate for use in lipid synthesis and energy generation.

The protein acts as a monomer and produces acetyl-CoA from acetate in a reaction that requires ATP.

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human, mouse

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu Mo

Immunogen:ACSS2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 38~68 amino acids from the N-terminal region of human ACSS2.

Shelf Life:6months

Swiss Prot Number:Q9NR19

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp