

US Biological (antibodypab)

84-1663-38 [取扱停止] GNAZ, NT (GNAZ, Guanine nucleotide-binding protein G(z) subunit alpha, G(x) alpha chain, Gz-alpha) (APC) 200ul 036136-APC



アズワン品番 84-1663-38

型番 036136-APC

入り数 1個

標準価格 -円

取扱停止



特徴

GNAZ is a member of a G protein subfamily that mediates signal transduction in pertussis toxin-insensitive systems.

This protein may play a role in maintaining the ionic balance of perilymphatic and endolymphatic cochlear fluids.

Applications: Suitable for use in Western Blot and ELISA.

Other applications not tested.

Recommended Dilutions: Optimal dilutions to be determined by the researcher.

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A and peptide affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Recognizes human GNAZ.

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu

Immunogen:KLH-conjugated synthetic peptide mapping to a fragment of residues with in amino acids 77-106 in the N-terminal region of human GNAZ.

Shelf Life:6months

Swiss Prot Number:P19086

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp