

AXEL

US Biological (antibodypab)

84-1663-47 [取扱停止] CLPTM1, ID (CLPTM1, Cleft lip and palate transmembrane protein 1) (APC) 200ul
034011-APC

USBiological
Life Sciences

アズワン品番 84-1663-47

型番 034011-APC

入り数 1個

標準価格 -円

取扱停止

USBiological
Life Sciences

特徴

CLPTM1 may play a role in T-cell development (By similarity).

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Recognizes human CLPTM1. Species Crossreactivity: mouse

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu Mo

Immunogen:KLH-conjugated synthetic peptide mapping to a fragment of residues with in amino acids 198-227 in the central region of human CLPTM1, UniProt Accession # O96005. Species sequence homology: bovine, zebrafish.

Shelf Life:6months

Swiss Prot Number:O96005

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp