

US Biological (antibodypab)

84-1663-58 [取扱停止] RPS11, ID (RPS11, 40S ribosomal protein S11) (APC) 200ul 041231-APC



アズワン品番 84-1663-58

型番 041231-APC

入り数 1個

標準価格 -円

取扱停止



特徴

Ribosomes, the organelles that catalyze protein synthesis, consist of a small 40S subunit and a large 60S subunit.

Together these subunits are composed of 4 RNA species and approximately 80 structurally distinct proteins.

This gene encodes a ribosomal protein that is a component of the 40S subunit.

The protein belongs to the S17P family of ribosomal proteins.

It is located in the cytoplasm.

The gene product of the E.

coli ortholog (ribosomal protein S17) is thought to be involved in the recognition of termination codons.

This gene is co-transcribed with a small nucleolar RNA gene, which is located in its third intron.

As is typical for genes encoding ribosomal proteins, there are multiple processed pseudogenes of this gene dispersed through the genome.

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size: 200ul

Host: rabbit

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity: Human

Isotype: IgG

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity: Hu

Immunogen: RPS11 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 87~116 amino acids from the Central region of human RPS11.

Shelf Life: 6 months

Swiss Prot Number: P62280

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp