

US Biological (antibodypab)

84-1663-62 [取扱停止] NPW, CT (NPW, PPL8, PPNP W, Neuropeptide W, Preproprotein L8, Neuropeptide W-23, Neuropeptide W-30) (APC) 200ul 039105-A PC



アズワン品番	84-1663-62
型番	039105-APC
入り数	1個

標準価格	-円
------	----

取扱停止



特徴

NPW is processed into 23- and 30-amino acid neuropeptides that bind and activate two G-protein coupled receptors in the central nervous system.

The neuropeptides have been shown to enhance cortisol secretion from adrenal cells through the adenylate cyclase/protein kinase A signaling cascade.

The preproprotein is translated using a non-AUG initiation codon that is inferred from analyses of the mouse ortholog.

[provided by RefSeq].

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu

Immunogen:NPW antibody is generated from rabbits immunized with a KLH conjugate d synthetic peptide between 128~158 amino acids from the C-terminal region of human NPW.

Shelf Life:6months

Swiss Prot Number:Q8N729

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp