



US Biological (antibodypab)

84-1663-72 [取扱停止] ERF3B, ID (GSPT2, ERF3B, Eukaryotic peptide chain release factor GTP-binding subunit ERF3B, G1 to S phase transition protein 2 homolog) (APC) 200ul 035201-APC



アズワン品番 84-1663-72

型番 035201-APC

入り数 1個

標準価格 -円

取扱停止



特徴

ERF3B encodes a GTPase that belongs to the GTP-binding elongation factor family. The encoded protein is a polypeptide release factor that complexes with eukaryotic peptide chain release factor 1 to mediate translation termination. This protein may also be involved in mRNA stability. Applications: Suitable for use in Western Blot, FLISA Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500 Storage and Stability: Store product at 4°C in the dark. DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid. Dilute required amount only prior to immediate use. Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human, mouse

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu Mo

Immunogen:ERF3B antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 223~253 amino acids from the Central region of human ERF3B.

Shelf Life:6months

Swiss Prot Number:Q8IYD1

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp