

US Biological (antibodypab)

84-1663-87 [取扱停止] DGAT2, NT (DGAT2, Diacylglycerol O-acyltransferase 2, Acyl-CoA retinol O-fatty-acyltransferase, Diglyceride acyltransferase 2) (APC) 200ul 034598-APC



アズワン品番 84-1663-87

型番 034598-APC

入り数 1個

標準価格 -円

取扱停止



特徴

Acyl-CoA:diacylglycerol acyltransferase, or DGAT (EC 2.3.1.20), is responsible for the synthesis of triglycerides.

It catalyzes a reaction in which diacylglycerol is covalently joined to long chain fatty acyl-CoAs.

Applications: Suitable for use in Western Blot and FLISA.

Other applications have not been tested.

Recommended Dilution: Optimal dilutions to be determined by the researcher.

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A and peptide affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Recognizes human DGAT2.

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu

Immunogen:KLH-conjugated synthetic peptide mapping to a fragment of residues with in amino acids 1-30 in the N-terminal region of human DGAT2, UniProt Accession #Q96PD7.

Shelf Life:6months

Swiss Prot Number:Q96PD7

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp