

US Biological (antibodypab)

84-1663-88 [取扱停止] GNB5, NT (Gbeta5, Transducin beta chain 5, Guanine nucleotide-binding protein subunit beta-5, GNB5) (APC) 200ul 036140-APC



アズワン品番 84-1663-88

型番 036140-APC

入り数 1個

標準価格 -円

取扱停止



特徴

Heterotrimeric guanine nucleotide-binding proteins (G proteins), which integrate signals between receptors and effector proteins, are composed of an alpha, a beta, and a gamma subunit.

These subunits are encoded by families of related genes.

This gene encodes a beta subunit.

Beta subunits are important regulators of alpha subunits, as well as of certain signal transduction receptors and effectors.

Applications: Suitable for use in Western Blot, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human, mouse

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu Mo

Immunogen:GNB5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 10~39 amino acids from the N-terminal region of human GNB5.

Shelf Life:6months

Swiss Prot Number:O14775

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp