



US Biological (antibodypab)

84-1663-94 [取扱停止] INSC, CT (INSC, Protein inscuteable homolog) (APC) 200ul 037056-APC



アズワン品番 84-1663-94

型番 037056-APC

入り数 1個

標準価格 -円

取扱停止



特徴

In Drosophila, neuroblasts divide asymmetrically into another neuroblast at the apical side and a smaller ganglion mother cell on the basal side.

Cell polarization is precisely regulated by 2 apically localized multiprotein signaling complexes that are tethered by Inscuteable, which regulates their apical localization (Izaki et al., 2006 [PubMed 16458856]).

Applications: Suitable for use in Western Blot, Immunohistochemistry, FLISA

Recommended Dilution: FLISA: 1:1,000 Western Blot: 1:100-500 Immunohistochemistry: 1:50-100

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Human, mouse

Isotype:IgG

Calc Applications Abbrev:FLISA IHC WB

Calc Crossreactivity:Hu Mo

Immunogen:INSC antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 416-445 amino acids from the C-terminal region of human INSC.

Shelf Life:6months

Swiss Prot Number:Q1MX18

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp