

US Biological (antibodypab)

84-1663-97 GDA, NT (GDA, KIAA1258, Guanine deaminase, Guanine aminohydrolase, p51-nedasin) (APC) 200ul 035976-APC



アズワン品番	84-1663-97
型番	035976-APC
入り数	1個
標準価格	-円



特徴

GDA is an enzyme that catalyzes the hydrolytic deamination of guanine, producing xanthine and ammonia.

Applications: Suitable for use in Western Blot, Flow Cytometry (Not Tested) and FLISA.

Other applications have not been tested.

Recommended Dilutions: Optimal dilutions to be determined by the researcher.

Storage and Stability: Store product at 4°C in the dark.

DO NOT FREEZE! Stable at 4°C for 12 months after receipt as an undiluted liquid.

Dilute required amount only prior to immediate use.

Further dilutions can be made in assay buffer.

Caution: APC conjugates are sensitive to light.

For maximum recovery of product, centrifuge the original vial prior to removing the cap.

Note: Applications are based on unconjugated antibody.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by Protein A and peptide affinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.2, 0.09% sodium azide and a stabilizer. Labeled with Allophycocyanin (APC).

Specificity:Recognizes human GDA.

Isotype:IgG

Calc Applications Abbrev:FLISA WB

Calc Crossreactivity:Hu

Immunogen:KLH-conjugated synthetic peptide mapping to a fragment of residues with in amino acids 94-122 in the N-terminal region of human GDA.

Shelf Life:6months

Swiss Prot Number:Q9Y2T3

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp