



US Biological (antibodypab)

84-4635-45 [取扱停止] RBBP8 200ul 393753



アズワン品番 84-4635-45

型番 393753

入り数 1個

標準価格 -円

取扱停止



特徴

The protein encoded by this gene is a ubiquitously expressed nuclear protein.

It is found among several proteins that bind directly to retinoblastoma protein, which regulates cell proliferation.

This protein complexes with transcriptional co-repressor CTBP.

It is also associated with BRCA1 and is thought to modulate the functions of BRCA1 in transcriptional regulation, DNA repair, and/or cell cycle checkpoint control.

It is suggested that this gene may itself be a tumor suppressor acting in the same pathway as BRCA1.

Three transcript variants encoding two different isoforms have been found for this gene.

More transcript variants exist, but their full-length natures have not been determined.

Applications: Suitable for use in ELISA, Immunohistochemistry.

Other applications not tested.

Recommended Dilution: Immunohistochemistry: 1:100-1:300 Optimal dilutions to be determined by the researcher.

Storage and Stability: May be stored at 4°C for short-term only.

Aliquot to avoid repeated freezing and thawing.

Store at -20°C.

Aliquots are stable for 12 months after receipt.

For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

仕様

Size:200ul

Host:rabbit

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by immunoaffinity chromatography.

Form:Supplied as a liquid in PBS, pH 7.3, 0.05% sodium azide, 50% glycerol.

Specificity:Recognizes human RBBP8. Species Crossreactivity: mouse and rat

Isotype:IgG

Calc Applications Abbrev:E IHC

Calc Crossreactivity:Hu Mo Rt

Immunogen:Recombinant protein corresponding to human RBBP8.

Shelf Life:1year

Swiss Prot Number:Q99708

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp