

AXEL

Abbexa

Kinesin-Like Protein KIF2C Phospho-Ser95 (KIF2C pS95) Antibody 50 ug

abbexa 

シリーズ共通仕様

Clonality	Polyclonal	Tested Applications	ELISA, WB, IF/ICC
Host	Rabbit	Conjugation	Unconjugated
Reactivity	Human, Mouse		

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
--------	--------	-----	----	----	------------------------------	----

 **abbexa** 

84-9237-07

Kinesin-Like Protein KIF2C Phospho-Ser95 (KIF2C pS95) Antibody 50 ug

abx325669

1個

75,000円

 **abbexa** 

84-9237-08

Kinesin-Like Protein KIF2C Phospho-Ser95 (KIF2C pS95) Antibody 100 ug

abx325669

1個

97,000円