




NeoBiotechnologies

CD20/ MS4A1 (B-Cell Marker) 100ug

特徴







Formulation:Purified Ab with BSA and Azide at 200ug/ml



シリーズ共通仕様

免疫動物	Mouse	交差性	Human.
別名	APY; ATOPY; B-lymphocyte cell-surface antigen B1; Bp35; Fc epsilon receptor I beta chain; Fc Fragment of IgE high affinity I receptor for beta polypeptide; FCER1B; High affinity immunoglobulin epsilon receptor subunit beta; IgE Fc receptor subunit beta; IGEL ; IGER; IGHER; LEU16; Leukocyte surface antigen Leu-16; Ly44; Membrane spanning 4 domains subfamily A member 2; Membrane-spanning 4-domains subfamily A member 1 (MS4A1)	Entrez Gene ID	931

Swiss -Prot	P11836	Unige ne	71 25 53
Hum an Ge ne Sy mbol	MS4A1	Hu-ch romos ome l ocati on	1 1 q 12. 2
Mol. weigh t of antig en	33-37kDa		

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	85-1237-66	CD20/ MS4A1 (B-Cell Marker) 20ug	931-MSM11-P0	1個	0円	
	85-1237-67	CD20/ MS4A1 (B-Cell Marker) 100ug	931-MSM11-P1	1個	0円	
	85-1237-68	CD20/ MS4A1 (B-Cell Marker) 100ug	931-MSM11-P1ABX	1個	0円	
	85-1294-09	CD20/ MS4A1 (B-Cell Marker) 20ug	931-MSM9-P0	1個	0円	
	85-1294-10	CD20/ MS4A1 (B-Cell Marker) 100ug	931-MSM9-P1	1個	0円	
	85-1294-11	CD20/ MS4A1 (B-Cell Marker) 100ug	931-MSM9-P1ABX	1個	0円	