

US Biological (nonantibody)

85-1663-77 Hydroxy Bosentan-d4 1mg CAS No:106

5472-91-4 013518

USBiological
Life Sciences

アズワン品番 85-1663-77

型番 013518

入り数 1個

標準価格 -円

USBiological
Life Sciences

特徴

- A labeled analogue of the active metabolite of Bosentan, a mixed endothelin receptor antagonist which is used as a vasodilator.
- Form:White Solid
- Melting Point:124-127°C
- Solubility:Methanol

仕様

- Size: 1mg
- Source Antigen: Synthetic
- Grade: Highly Purified
- Form: White Solid
- Molecular Formula: C₂₇H₂₅D₄N₅O₇S
- Molecular Weight Non Ab: 571.64
- EU Commodity Code: 38220090
- この商品は法規制を確認していません。(法規制によって販売できない場合もございます)
- 製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。
- CAS No : 1065472-91-4

- 【試薬に関するお問合せ】
- アズワン株式会社 試薬・プロセス材料グループ
- TEL: 06-6447-8641
- FAX: 06-6447-8642
- E-mail: Reagents@so.as-1.co.jp