

AXEL

US Biological (nonantibody)

85-1663-77 Hydroxy Bosentan-d4 1mg 013518

USBiological
Life Sciences

アズワン品番 85-1663-77

型番 013518

入り数 1個

標準価格 -円

USBiological
Life Sciences

特徴

A labeled analogue of the active metabolite of Bosentan, a mixed endothelin receptor antagonist which is used as a vasodilator.

Form:White Solid

Melting Point:124-127°C

Solubility:Methanol

仕様

Size:1mg

Source Antigen:Synthetic

Grade:Highly Purified

Form:White Solid

Molecular Formula:C27H25D4N5O7S

Molecular Weight Non Ab:571.64

EU Commodity Code:38220090

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp