

US Biological (nonantibody)

85-1663-88 rac 3-Hydroxybutyric Acid-d4 Sodium Salt 5mg CAS No:1219804-68-8, Unlabeled salt: 150-83-4 013533

US Biological
Life Sciences

アズワン品番 85-1663-88

型番 013533

入り数 1個

標準価格 -円

US Biological
Life Sciences

特徴

- Labeled rac 3-Hydroxybutyric Acid. Optically active 3-hydroxybutyric acids are key intermediates of the biosynthesis and metabolism of fatty acids and exist widely in biological systems.
- Synonyms:3-Hydroxybutanoic Acid-d4 Sodium Salt; 3-Hydroxybutyric Acid-d4 Monosodium Salt; DL-β-Hydroxybutyric Acid-d4 Sodium Salt; Sodium 3-Hydroxybutyrate-d4; Sodium DL-β-Hydroxybutyrate-d4; Sodium β-Hydroxybutyrate-d4; β-Hydroxybutyric Acid-d4 Sodium Salt
- Molecular Formula:C₄H₇O₂aO₂
- Molecular Weight:130.11
- Purity:≥96%
- Appearance:White to off-white solid
- Melting Point:163-165°C (lit)
- Solubility:DMSO (Slightly, Heated), Methanol (Slightly)
- Method for Determining Identity:Proton NMR and Mass Spectrometric Analysis
- Elemental Analysis:As Reported
- Sodium Content (Ion Chromatography):As Reported
- Storage and Stability:May be stored at RT for short-term only. Long-term storage is recommended at -20°C for 6 months after receipt. Hygroscopic.

仕様

- Size:5mg
- Source Antigen:Synthetic
- Grade:Highly Purified
- Molecular Formula:C₄H₇O₂aO₂
- Molecular Weight Non Ab:130.11

- EU Commodity Code:38220090
- この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)
- 製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。
- CAS No :1219804-68-8,Unlabeled salt: 150-83-4

- 【試薬に関するお問合せ】
- アズワン株式会社 試薬・プロセス材料グループ
- TEL:06-6447-8641
- FAX:06-6447-8642
- E-mail:Reagents@so.as-1.co.jp