



US Biological (nonantibody)

Signal Transducer And Activator Of Transcription 6 (STAT6) Recombinant, Mouse 50ug



特徴

Source: Recombinant Mouse from E. coli

Purity: $\geq 95\%$

Endotoxin: 1.0 EU per 1 μ g (determined by the LAL method)

Accession No: P52633

Fragment: Leu410~Glu623 (Accession No: P52633)

Sequence: MGSSHHHHHH SSGLVPRGSH MASMTGGQ

QM GRGSEF-L VVIVHGNQDN NAKATILWDN AFSEMDR

VVF VVAERPWEK MCETLNLFKFM AEVGTSRGLL PEHFLF

LAQK IFNDNSLSVE AFQHRCVSWS QFNKEILLGR GFTFW

QWFDG VLDLTKRCLR SYWSDRLIIG FISKQYVTSL LLNE

PDGTFL LRFSDSEIGG ITIAHVIRGQ DGSSQIENIQ PFSAK

DLSIR SLGDRIRDLA QLKNLYPKKP KDE

Epitope Tag: N-terminal Tags: His-tag and T7-tag

Molecular Weight: 28.4 kD

Applications: Suitable for use in ELISA, Western Blot, Immunoprecipitation, SDS-PAGE

Other applications not tested.

Recommended Dilution: Optimal dilutions to be determined by the researcher.

Storage and Stability: May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -80°C. Aliquots are stable for at least 12 months from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

Note: Thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation by incubating the protein at 37°C for 48h, with no obvious degradation or precipitation.

ion observed. Protein loss is <5% within the expiration d
ate under appropriate storage conditions.

シリーズ共通仕様

Source Antigen	Recombi nant Mous e from E. c oli	Grade	Highly Purified
Purity	≥95%	Form	Supplied as lyophilized powder in PBS, pH7.4, 5% sucrose, 0. 01% sarcosyl. Reconstitute in s terile PBS, pH7.2.
Molecul ar Weight Non Ab	28.4	Swiss P rot N umber	P52633
EU Co mmodity Code	30021019		

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	85-1908-75	Signal Transducer And Activator Of Transcription 6 (STAT6) Recombina nt, Mouse 10ug	156815	1個	98,000円	
	85-1908-76	Signal Transducer And Activator Of Transcription 6 (STAT6) Recombina nt, Mouse 50ug	156815	1個	155,000 円	
	85-1908-77	Signal Transducer And Activator Of Transcription 6 (STAT6) Recombina nt, Mouse 200ug	156815	1個	232,000 円	