



US Biological (nonantibody)

Signal Transducer And Activator Of Transcription 6 (STAT6) Recombinant, R at 200ug



特徴

Source: Recombinant Rat from E. coli

Purity: $\geq 95\%$

Endotoxin: 1.0 EU per 1 μ g (determined by the LAL method)

Accession No: Q1KQ07

Fragment: Gly557~Trp841 (Accession No: Q1KQ07)

Sequence: MGSSHHHHHH SSGLVPRGSH MASMTGGQ
QM GRGSEF- GTFL LRFSDSEIGG ITIAHVIRGQ DGSSQIE
NIQ PFSAKDLSIR SLGDRIRDLA QLKNLYPKKP KDEAFRSH
YK PEQMGKDGRG YVSTTIKMTV ERDQPLPTPE PQMPA
MVAPY DLGMAPDASM QLSSDMVYPP QSHSIHSFQS IPLE
ESMSVL PPFQEPHLQM PPTMSQISMP FDQPHPQGGLL Q
CQSQEHA VS SPELLCSDV TMAEDSCLTQ PVQGFPQGT
W VSEGMYPPLM PPTEQDLTKL LLEGQGEGGG SIGTQPLL
QP SSYGQSGISM SHLDLRTNPS W

Epitope Tag: N-terminal Tags: His-tag and T7-tag

Molecular Weight: 35.1 kD

Applications: Suitable for use in ELISA, Western Blot, Immunoprecipitation, SDS-PAGE

Other applications not tested.

Recommended Dilution: Optimal dilutions to be determined by the researcher.

Storage and Stability: May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -80°C. Aliquots are stable for at least 12 months from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

Note: Thermal stability is described by the loss rate of the target protein. The loss rate was determined by acce

lerated thermal degradation by incubating the protein at
 37°C for 48h, with no obvious degradation or precipitat
 ion observed. Protein loss is <5% within the expiration d
 ate under appropriate storage conditions.

シリーズ共通仕様

Source Antigen	Recombi nant Rat from E. c oli	Grade	Highly Purified
Purity	≥95%	Form	Supplied as lyophilized powder in PBS, pH7.4, 5% sucrose, 0.01 % sarcosyl. Reconstitute in steri le PBS, pH7.2.
Molecul ar Weight Non Ab	35.1	Swiss P rot N umber	Q1KQ07
EU Co mmodity Code	3002101 9		

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	85-1908-78	Signal Transducer And Activator Of Transcription 6 (STAT6) Recombina nt, Rat 10ug	156816	1個	100,000 円	
	85-1908-79	Signal Transducer And Activator Of Transcription 6 (STAT6) Recombina nt, Rat 50ug	156816	1個	158,000 円	
	85-1908-80	Signal Transducer And Activator Of Transcription 6 (STAT6) Recombina nt, Rat 200ug	156816	1個	244,000 円	