



US Biological (nonantibody)

Signal Transducing Adaptor Protein 1 (STAP1) Recombinant, Mouse 200ug



特徴

Source: Recombinant Mouse from E. coli

Purity: $\geq 95\%$

Endotoxin: 1.0 EU per 1 ug (determined by the LAL method)

Accession No: Q9JM90

Fragment: Met1 ~ Ala297 (Accession No: Q9JM90)

Sequence: MGSSHHHHHH SSGLVPRGSH MASMTGGQ
QM GRGS- MMAKKPPKPA PRRIFQERLK ITALPLYFEG FLL
VKRSDHQ EYKHYWTELR GTTLFFYTDK KSTIYVGKLD II
DLVCLTGQ HSTEKNCAKF TLVLPKEEVH VKTENTESGE E
WRGFILTVT ELTVPQHVS L LPGQVIRLHE VLEREKKRRI ET
DQLPLMPP EKEKEPVQDY ADVLNPLPEC FYAVSRKEAT A
MLEKNPSWG NMILRPGSDS KNYSITIRQE IEMPRIKHFK
VTRTGNNYTI ELEKPVTLPN LFSVIDYFVK ETRGNLRPFI H
SADDNFGQD PNIEDRSEKF KKNPHNA

Epitope Tag: N-terminal Tags: His-tag and T7-tag

Molecular Weight: 38.2 kD

Applications: Suitable for use in ELISA, Western Blot, Immunoprecipitation, SDS-PAGE

Other applications not tested.

Recommended Dilution: Optimal dilutions to be determined by the researcher.

Storage and Stability: May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -80°C. Aliquots are stable for at least 12 months from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

Note: Thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation by incubating the protein at

37°C for 48h, with no obvious degradation or precipitation observed. Protein loss is <5% within the expiration date under appropriate storage conditions.

シリーズ共通仕様

Source Antigen	Recombinant Mouse from E. coli	Grade	Highly Purified
Purity	≥95%	Form	Supplied as lyophilized powder in PBS, pH7.4, 5% sucrose, 0.01% sarcosyl. Reconstitute in sterile PBS, pH7.2.
Molecular Weight Non Ab	38.2	Swiss Protein Number	Q9JMQ0
EU Commodity Code	30021019		

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	85-1908-81	Signal Transducing Adaptor Protein 1 (STAP1) Recombinant, Mouse 10 ug	156817	1個	98,000円	
	85-1908-82	Signal Transducing Adaptor Protein 1 (STAP1) Recombinant, Mouse 50 ug	156817	1個	155,000 円	
	85-1908-83	Signal Transducing Adaptor Protein 1 (STAP1) Recombinant, Mouse 20 0ug	156817	1個	232,000 円	