



US Biological (nonantibody)

Signaling Threshold Regulating Transmembrane Adaptor 1 (SIT1) Recombinant, Rat 50ug



特徴

Source: Recombinant Rat from E. coli

Purity: $\geq 95\%$

Endotoxin: 1.0 EU per 1 μ g (determined by the LAL method)

Accession No: Q5M869

Fragment: Ser2~Pro175 (Accession No: Q5M869)

Sequence: MGSSHHHHHH SSGLVPRGSH MASMTGGQ
QM GRGSEF- SRENNCTTA DLAWGIPSIT QAWGLWALFG
VVTMLLLISL AALLSQWTRG RRRTQEEQGP PSGRSVEEVP
LYGNLHYLQT GRLSEESRSE EQDPSSGGLA RGAEAEATCYT
SLQLRPAQGR IPSSGTPIKY CEVVDSEPK PQASGPEPEL Y
ASVCAQTRR ARASFPDQAY ANSQP

Epitope Tag: N-terminal Tags: His-tag and T7-tag

Molecular Weight: 22.8 kD

Applications: Suitable for use in ELISA, Western Blot, Immunoprecipitation, SDS-PAGE

Other applications not tested.

Recommended Dilution: Optimal dilutions to be determined by the researcher.

Storage and Stability: May be stored at 4°C for short-term only. Aliquot to avoid repeated freezing and thawing. Store at -80°C. Aliquots are stable for at least 12 months from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

Note: Thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation by incubating the protein at 37°C for 48h, with no obvious degradation or precipitation observed. Protein loss is <5% within the expiration date.

ate under appropriate storage conditions.

シリーズ共通仕様

Source Antigen	Recombi nant Rat from E. c oli	Grade	Highly Purified
Purity	≥95%	Form	Supplied as lyophilized powder in PBS, pH7.4, 5% sucrose, 0.01 % sarcosyl. Reconstitute in steri le PBS, pH7.2.
Molecul ar Weight Non Ab	22.8	Swiss P rot N umber	Q5M869
EU Co mmodity Code	3002101 9		

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
	85-1908-90	Signaling Threshold Regulating Tr ansmembrane Adaptor 1 (SIT1) Re combinant, Rat 10ug	156820	1個	99,000円	
	85-1908-91	Signaling Threshold Regulating Tr ansmembrane Adaptor 1 (SIT1) Re combinant, Rat 50ug	156820	1個	157,000 円	
	85-1908-92	Signaling Threshold Regulating Tr ansmembrane Adaptor 1 (SIT1) Re combinant, Rat 200ug	156820	1個	241,000 円	