

AXEL

US Biological (nonantibody)

85-4635-45 Mothers Against Decapentaplegic Homolog 2 (Smad2) BioAssay™ ELISA Kit (Human) 96Tests 366328



アズワン品番	85-4635-45
型番	366328
入り数	1個
標準価格	-円



特徴

Sample Type:Tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Intended Use:BioAssay™ ELISA Kit for Mothers Against Decapentaplegic Homolog 2 (Smad2) in tissue homogenates, cell lysates, cell culture supernates and other biological fluids

Range:0.156-10ng/ml

Specificity:Human

Storage and Stability:Store powder at 4°C liquid at -20°C. Store other components at 4°C. Stable for at least 6 months For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

仕様

Size:96Tests

Source Antigen:Human

Calc Applications Abbrev:E

EU Commodity Code:38220000

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp