



US Biological (nonantibody)

85-4635-87 Nicotinic Acid Adenine Dinucleotide Phosphate (NAADP) BioAssay™ ELISA Kit (General) 96 Tests 366370



アズワン品番 85-4635-87

型番 366370

入り数 1個

標準価格 -円



特徴

Sample Type: Serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.

Intended Use: BioAssay™ ELISA Kit for Nicotinic Acid Adenine Dinucleotide Phosphate (NAADP) in serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.

Range: 123.46-10,000ng/ml

Specificity: General

Storage and Stability: Store powder at 4°C liquid at -20°C. Store other components at 4°C. Stable for at least 6 months For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap.

仕様

Size:96Tests

Source Antigen:General

Calc Applications Abbrev:E

EU Commodity Code:38220000

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp