



US Biological (nonantibody)

4-chloro-3,5-dimethoxy-3'-quinolinecarbonitrile 50mg

特徴

4-chloro-3,5-dimethoxy-3'-quinolinecarbonitrile is used as a reactant in the synthesis of Src kinase activity inhibitors.



Molecular Formula:C12H9ClN2O2

Molecular Weight:248.67

Storage and Stability:Store at -20°C. For maximum recovery of product, centrifuge the original vial prior to removing the cap.

シリーズ共通仕様

Grade	Highly Purified	Molecular Formula	C12H9ClN2O2
Molecular Weight Non Ab	248.67	EU Commodity Code	38220090

商品イメージ	アズワン品番	商品名	型番	入数	標準価格 (税抜) 定額送料 (税抜)	備考
--------	--------	-----	----	----	------------------------------	----



85-5259-83

4-chloro-3,5-dimethoxy-3'-quinolinecarbonitrile 50mg

431777

1個

97,000円



85-5259-84

4-chloro-3,5-dimethoxy-3'-quinolinecarbonitrile 250mg

431777

1個

153,000円