

AXEL

US Biological (antibodymab)

85-8317-07 Interleukin 1 Family Member 7 (IL-1F7, IL1F7, Interleukin-1 zeta, IL-1 zeta, FIL1 zeta, FIL1(zeta), FIL1Z, FIL1, Interleukin-1 Homolog 4, IL-1H4, IL1H4, Interleukin-1-related Protein, IL-1RP1, IL1RP1, IL-1X, Interleukin 1 Family, Member 7 (zeta)) (Azide Free) (MaxLight 490) 100ul I7663-24R6-ML490

USBiological
Life Sciences

アズワン品番	85-8317-07
型番	I7663-24R6-ML490
入り数	1個
標準価格	-円

USBiological
Life Sciences

仕様

Size:100ul

Host:mouse

Source Antibody:human

Grade:Affinity Purified

Purity:Purified by immunoaffinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

490.

Specificity: Recognizes human IL1F7.

Clone # USB: 11C309

Isotype: IgG1

Calc Applications Abbrev: IHC

Calc Crossreactivity: Hu

Immunogen: Full length recombinant protein corresponding to aa1-219 of human IL1F7, with GST tag. MW of GST tag alone is 26kD.

Shelf Life: 1 year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp