

AXEL

US Biological (antibodymab)

85-8317-29 SPAG6 (Sperm-Associated Antigen 6, Protein PF16 Homolog, Sperm Flagellar Protein, Repro-
-SA-1, PF16) (Azide Free) (MaxLight 490) 100ul S234
3-ML490

USBiological
Life Sciences

アズワン品番	85-8317-29
型番	S2343-ML490
入り数	1個
標準価格	-円

USBiological
Life Sciences

仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Purified

Purity: Purified

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 490.

Specificity: Recognizes human SPAG6 at ~55kD.

Clone # USB: 11C774

Isotype:IgG3,k

Calc Applications Abbrev:WB

Calc Crossreactivity:Hu

Immunogen:Recombinant protein corresponding to a portion of amino acids within aa 1-101 of human SPAG6 with a GST Tag (MW of tag alone is 26kD). (O75602, NP_036575)

Shelf Life:1year

Swiss Prot Number:O75602

EU Commodity Code:30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp