



US Biological (antibodymab)

85-8317-34 [取扱停止] Tankyrase (TRF1-interacting Ankyrin-related ADP-ribose Polymerase, TNKS, Tankyrase I, Tankyrase-I, TANK1, TNKS1, TNKS-1, PARP-5a, PARP5a, PARPL, TIN1, TINF1) (Azide Free) (MaxLight 490) 100ul T0900-20B-ML490



アズワン品番	85-8317-34
型番	T0900-20B-ML490
入り数	1個
標準価格	-円
取扱停止	



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified from ascites by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

tTM490.

Specificity: Recognizes human TNKS.

Clone # USB: 11C843

Isotype: IgG2a

Calc Applications Abbrev: FLISA IHC WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant protein corresponding to a region near a portion of amino acids within aa1017-1125 of human TNKS with a GST tag. (NP_003738)

Shelf Life: 1year

Swiss Prot Number: O95271

EU Commodity Code: 30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp