

US Biological (antibodymab)

85-8317-45 WRB (Tryptophan-rich Protein, Congenital Heart Disease 5 Protein, CHD5, WRB Tryptophan Rich Basic Protein) (Azide Free) (MaxLight 490) 100 ul W8900-ML490



アズワン品番	85-8317-45
型番	W8900-ML490
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Purified

Purity: Purified

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 490.

Specificity: Recognizes human WRB.

Clone # USB: 11C959

Isotype: IgG2a,k

Calc Applications Abbrev:WB

Calc Crossreactivity:Hu

Immunogen:Partial recombinant protein corresponding to a portion of amino acids within aa29-102 of human WRB (NP_004618) with GST tag. MW of GST tag alone is 26kD.

Shelf Life:1year

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp