



US Biological (antibodymab)

85-8317-83 [取扱停止] ADA, CT (ADA1, Adenosine deaminase, Adenosine aminohydrolase) (MaxLight 490) 100ul 030746-ML490



アズワン品番	85-8317-83
型番	030746-ML490
入り数	1個
標準価格	-円

取扱停止



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein G affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 490.

Specificity: Recognizes human ADA.

Clone # USB: 608CT2.1.3

Isotype: IgG1

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity:Hu

Immunogen:KLH-conjugated synthetic peptide mapping to a fragment of residues with
in amino acids 287-314 in the C-terminal region of human ADA.

Shelf Life:1year

Swiss Prot Number:P00813

EU Commodity Code:30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL:06-6447-8641

FAX:06-6447-8642

E-mail:Reagents@so.as-1.co.jp