

AXEL

US Biological (antibodymab)

85-8317-87 [取扱停止] AKT3 (AKT3, PKBG, RAC-gamma serine/threonine-protein kinase, Protein kinase Akt-3, Protein kinase B gamma, RAC-PK-gamma, STK-2) (MaxLight 490) 100ul 030758-ML490

US Biological
Life Sciences

アズワン品番	85-8317-87
型番	030758-ML490
入り数	1個

標準価格	-円
------	----

取扱停止

US Biological
Life Sciences

仕様

Size: 100ul

Host: rabbit

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein G affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 490.

Specificity: Recognizes human AKT3. Species Crossreactivity: mouse

Clone # USB: 154CT4.5.3

Isotype: IgG

Calc Applications Abbrev: FLISA IHC WB

Calc Crossreactivity: Hu Mo

Immunogen: Recombinant protein corresponding to mouse immunized with AKT3. Uni

Prot Accession: Q9Y243

Shelf Life: 1year

Swiss Prot Number: Q9Y243

EU Commodity Code: 30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp