

US Biological (antibodymab)

85-8564-13 PHKB (Phosphorylase B Kinase Regulatory Subunit beta, Phosphorylase Kinase Subunit beta, DKFZp781E15103, FLJ41698) (MaxLight 650) 100ul 131268-ML650



アズワン品番	85-8564-13
型番	131268-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight 650.

Specificity: Recognizes human PHKB.

Clone # USB: 2E9

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa984-1094 from human PHKB (NP_000284) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp