

US Biological (antibodymab)

85-8564-14 PHKG1 (Phosphorylase Kinase Subunit
gamma 1, Phosphorylase B Kinase gamma Catalytic C
hain Skeletal Muscle Isoform, PHKG) (MaxLight 650)
100ul 131269-ML650



アズワン品番	85-8564-14
型番	131269-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLigh

tTM650.

Specificity: Recognizes human PHKG1.

Clone # USB: 2B4

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA IHC WB

Calc Crossreactivity: Hu

Immunogen: Full length recombinant corresponding to aa1-388 from human PHKG1 (AAH51327) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認しておりません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp