

US Biological (antibodymab)

85-8564-16 PHLDA2 (Pleckstrin Homology-like Domain Family A Member 2, Beckwith-Wiedemann Syndrome Chromosomal Region 1 Candidate Gene C Protein, Imprinted in Placenta and Liver Protein, Tumor-suppressing STF cDNA 3 Protein, Tumor-suppressing Subchromosomal Transferable Fragment Candidate Gene 3 Protein, p17-Beckwith-Wiedemann Region 1 C, p17-BWR1C, BWR1C, HLDA2, IPL, TSSC3) (MaxLight 650) 100ul 131276-ML650

USBiological
Life Sciences

アズワン品番	85-8564-16
型番	131276-ML650
入り数	1個
標準価格	-円

USBiological
Life Sciences

仕様

Size:100ul

Host:mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight™650.

Specificity: Recognizes human PHLDA2.

Clone # USB: 5E3

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa1-110 from human PHLDA2 (NP_003302) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp