

US Biological (antibodymab)

85-8564-20 Phosphodiesterase 9A (PDE9A, High Affinity cGMP-specific 3',5'-cyclic Phosphodiesterase 9A, FLJ90181, HSPDE9A2) (MaxLight 650) 100ul 131282-ML650



アズワン品番	85-8564-20
型番	131282-ML650
入り数	1個
標準価格	-円



仕様

Size: 100ul

Host: mouse

Source Antibody: human

Grade: Affinity Purified

Purity: Purified by Protein A affinity chromatography.

Form: Supplied as a liquid in PBS, pH 7.2. No preservative added. Labeled with MaxLight

650.

Specificity: Recognizes human PDE9A.

Clone # USB: 1E1

Isotype: IgG2a,k

Calc Applications Abbrev: FLISA IHC WB

Calc Crossreactivity: Hu

Immunogen: Partial recombinant corresponding to aa434-534 from human PDE9A (NP_001001567) with GST tag. MW of the GST tag alone is 26kD.

Shelf Life: 1year

EU Commodity Code: 30021010

この商品は法規制を確認していません。(法規制によって販売できない場合もございます)

製品の仕様は予告なく変更になる場合がございます。最新仕様はメーカーホームページをご確認ください。

【試薬に関するお問合せ】

アズワン株式会社 試薬・プロセス材料グループ

TEL: 06-6447-8641

FAX: 06-6447-8642

E-mail: Reagents@so.as-1.co.jp